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## LITERATUR

Al Zahran K, Al Jahdali H, Poirier L, Rene P, Gennaro ML, Menzies D. Accuracy and utility of commercially available amplification and serologic tests for the diagnosis of minimal pulmonary tuberculosis. *Am J Respir Crit Care Med.* 2000 Oct; 162 (4 Pt 1): 1323-9.

Baumann R, Kaempfer S, Chegou NN, Nene NF, Veenstra H, Spallek R, Bolliger CT, Luky PT, van Heiden PD, Singh M, Walzl G. Serodiagnostic markers for the prediction of the outcome of intensive phase tuberculosis therapy. *Tuberculosis (Edinb).* 2013 Mar;93(2):239-45. doi: 10.1016/j.tube.2012.09.003. Epub 2012 Nov 2.

Baumann R, Kaempfer S, Chegou N N., Oehlmann W, Spallek R, Loxton A G.; van Heiden P D.; Black, Gillian F.; Singh, M.; Walzl, G. 2015. A Subgroup of Latently Mycobacterium tuberculosis Infected Individuals Is Characterized by Consistently Elevated IgA Responses to Several Mycobacterial Antigens. *Mediators of Inflammation.* Vol. 2015, p1-10. 10p

Delphj, Chatterjee and Kay-Hooi Khoo: Mycobacterial lipoarabinomannan: An extraordinary lipoheteroglycan with profound physiological effects. *Glycobiology.* 1998 Feb; 8(2):113-20.

Dolapo Awoniyi, Raff Baumann, Novel Chegou, Belinda Kriel, Ruschca Jacobs, Martin Kidd, Andre Loxton, Susanne Kaempfer, Mahavir Singh, Gerhard Walzl. Combined specific IgG – and IgA –based diagnosis of Tuberculosis in African primary healthcare clinic attendees with signs and symptoms suggestive of Tuberculosis. DOI: 10.1136/bmjgh-2016-000260.94 Published 12 February 2017

Friscia G, Vordermeier HM, Pasvol G, Harris DP, Moreno C, Ivanyi: Human T cell responses to peptide epitopes of the 16-kD antigen in tuberculosis. *J. Clin Exp Immunol.* 1995 Oct; 102 (1): 53-7. Verbon A, Gennaro ML. Immunologic diagnosis of tuberculosis. *Clin Infect Dis.* 2000 Jun; 30 Suppl 3: 243-6. Review.

Imaz MS, Comini MA, Zerbin E, Sequeira MD, Spoletti MJ, Etchart AA, Pagano HJ, Bonifasch E, Diaz N, Claus JD, Singh M. Evaluation of the diagnostic value of measuring IgG, IgM and IgA antibodies to the recombinant 16-kilodalton antigen of mycobacterium tuberculosis in childhood tuberculosis. *Int J Tuberc Lung Dis* 2001 5 (11): 1036-43

Sardella I G; Mahavir S; Kumpfer S; Ribeiro Heringer R ; Féres Saad M H, +; Puccioni Sohier M. Evaluation of Lionex TB kits and mycobacterial antigens for IgG and IgA detection in cerebrospinal fluid from tuberculosis meningitis patients. Print version ISSN 0074-0276, Mem.

Kaisermann MG, Sardella IG, Trajman A, Coelho LV, Kämpfer S, Jonas F, Singh M, Saad MH.: IgA antibody responses to Mycobacterium tuberculosis recombinant MPT-64 and MT-10.3 (Rv3019c) antigens in pleural fluid of patients with tuberculous pleurisy. *Int J Tuberc Lung Dis.* 2005 Apr;9(4):461-6. PMID: 15830754

Khoo, K. H., Tang, J. B., and Chatterjee, D.: Variation in Mannose-capped Terminal Arabinan Motifs of Lipoarabinomannans from Clinical Isolates of Mycobacterium tuberculosis and Mycobacterium avium Complex (2001) *J. Biol. Chem.* 276, 3863-3871

Kuijper S, Jansen HM, Speelman P, Kolk AH: Development of a serological test for tuberculosis *Ned Tijdschr Geneeskd* 1991 Jan 26; 135 (4): 134-8

Lenka M, Pereira Arias-Bouda, Lan N, Nguyen, Ly M, Ho, Sjoukje Kuijper, Henk M, Jansen, and Arend H. J. Kolk: Development of Antigen Detection Assay for Diagnosis of Tuberculosis Using Sputum Samples *Journal of Clinical Microbiology*, June 2000, p. 2278-2283, Vol.38, No.6 0095-1137/00/504.00-0

Lyschchenko KP, Singh M, Colangeli R, Gennaro ML. A multi-antigen print immunoassay for the development of serological diagnosis of infectious diseases. *J Immunol Methods.* 2000 Aug 28; 242 (1-2): 91-100.

MURAMIRA M. Nobert, Dr OLAF MÜLLER, Dr. MUSEMAKWERI André: EVALUATION OF TB-ST TUBERCULOSIS RAPID TEST IN ADULT PULMONARY TUBERCULOSIS IN RWANDA. Huye, Université Nationale du Rwanda, FACULTY OF MEDICINE October 2007.

Pant Pai, N and Pai M: Point-of-Care Diagnostics for HIV and Tuberculosis: Landscape, Pipeline, and unmet needs. *Discovery Medicine*, Volume 13, Number 68 p.35-45. Jan. 2012. ISSN: 1539-6509

Pukazhvanthen P, Anbarasu D, Ahmed Kabeer Basirudeen S, Raja A, \*, Singh M: Assessing humoral immune response of 4 recombinant antigens for serodiagnosis of tuberculosis. *Tuberculosis* 94 (2014) 622e633.

Samanich KM, Keen MA, Vissa VD, Harder JD, Spencer JS, Belisle JT, Zolla-Pazner S, Laal S. Serodiagnostic potential of culture filtrate antigens of Mycobacterium tuberculosis. *Clin Diagn Lab Immunol.* 2000 Jul; 7(4): 662-8.

Samanich K, Belisle JT, Laal S. Homogeneity of antibody responses in tuberculosis patients. *Infect Immun.* 2001 Jul; 69 (7): 4600-9.

Silva VM, Kanaujia G, Gennaro ML, Menzies D. Factors associated with humoral response to ESAT-6, 38 kDa and 14 kDa in patients with a spectrum of tuberculosis. *Int J Tuberc Lung Dis.* 2003 May; 7 (5): 478-84.

Van Deun A, Portaels, F: Limitations and requirements for quality control of sputum smear microscopy for acid-fast bacilli. *J Clin Microbiol* 2000 Jun;38(6):2278-83.

Verbon A, Kuijper S, Jansen H, Speelman P, Kolk, A H J. Antigens in culture supernatant of *M. tuberculosis*: epitopes defined by monoclonal and human antibodies. *J Gen Microbiol.* 1990;136:955-964

Wilkinson, R. J., K. Hasløv, R. Rappuoli, F. Giovannoni, P. R. Narayanan, C. R. Desai, M. Vordermeier, J. Paulsen, G. Pasvol, J. Ivanyi, and M. Singh. Evaluation of the recombinant 38-kilodalton antigen of Mycobacterium tuberculosis as a potential immunodiagnostic reagent. *J. Clin. Microbiol.* 1997. 35: 553-557

WHO Tuberculosis Diagnostics Workshop: – Product Development Guidelines – Cleveland, Ohio, 27 July, 1997.

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