COVID-19 Ag Rapid Test / COVID-19 Ag Schnelltest

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Σ 20 Tests

REF COV19_AG_20

The COVID-19 Ag Rapid Test is a visual test for the direct and qualitative detection of SARS-CoV-2 viral Spike glycoprotein (S1) antigen in human nasopharynx and oropharynx within 15 min (latest 20 min).

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
LIONEX Diagnostics
and Therapeutics

Content / Inhalt	
INTENDED USE	3
INTRODUCTION / FIELD OF APPLICATION	3
PRINCIPLE OF THE TEST	3
SUPPLIED MATERIALS	3
MATERIALS NEEDED BUT NOT SUPPLIED	3
PREPARATION OF REAGENTS	4
STABILITY AND STORAGE CONDITIONS	4
WARNINGS AND PRECAUTIONS	4
SAMPLE COLLECTION AND PREPARATION	5
TEST PROCEDURE	5
INTERPRETATION OF RESULTS	7
QUALITY CONTROL	7
PERFORMANCE CHARACTERISTICS	7
Interfering Substances	9
Cross reactivity	9
LIMITATIONS	10
VERWENDUNGSZWECK	13
ZUSAMMENFASSUNG / ANWENDUNGSBEREICH	13
TESTPRINZIP	13
MITGELIEFERTE MATERIALIEN	13
ZUSÄTZLICHE ERFORDERLICHE MATERIALIEN, NICHT IM LIEFERUMFANG ENTHALTEN	13
VORBEREITUNG VON REAGENZIEN	14
STABILITÄTS- UND LAGERBEDINGUNGEN	14
WARNUNGEN UND VORSICHTSMAßNAHMEN	14
PROBENENTNAHME UND -VORBEREITUNG	15
TESTDURCHFÜHRUNG	15
INTERPRETATION DER ERGEBNISSE	17
QUALITÄTSKONTROLLE	17
LEISTUNGSMERKMALE	17
Störende Substanzen	19
Kreuzreaktivität	19
EINSCHRÄNKUNGEN	20
LITERATURE / LITERATUR	23
SYMBOLS / SYMBOLE	24
CONTACT INFORMATION / KONTAKTINFORMATIONEN	24

INTENDED USE

The **COVID-19 Ag Rapid Test** is a visual test for the direct and qualitative detection of SARS-CoV-2 viral Spike glycoprotein (S1) antigen in human nasopharynx and oropharynx within 15 min (latest 20 min). The test is for professional in vitro diagnostic use only and intended to aid in the diagnosis of SARS-CoV-2 infection and to complement direct pathogen detection using molecular methods. The product is intended for use as an IVD but can also be used for research purposes.

INTRODUCTION / FIELD OF APPLICATION

In December 2019, a novel zoonotic coronavirus SARS-CoV-2 was identified as an infectious agent that could cause an outbreak of viral pneumonia in human. Common signs of infection with the coronavirus include respiratory symptoms, breathing difficulties, fever, sore throat, stuffy nose, and dry cough. In some severe cases, the infection can cause viral pneumonia, severe acute respiratory syndrome (SARS), as well kidney failure and finally death. To prevent an infection with the coronavirus, it is recommended to avoid close contact with anyone showing symptoms of respiratory illness. One should follow standard hygienic methods like hand washing, as well covering mouth and nose. SARS-CoV-2 has structure proteins including spike (S), envelope (E), membrane (M) and nucleocapsid (N). The spike protein (S) is a glycoprotein, composed of two subunits (S1 and S2). It was found that the subunit S1 contains a receptor binding domain (RBD) which strongly interacts with human ACE2 receptor, causing infection of human respiratory cells.

COVID-19 Ag Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2 S1 from swabs. This test is for in-vitro professional diagnostic use and intended as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms and signs consistent with COVID-19. It provides only an initial screening test result and more specific alternative diagnosis methods should be performed to obtain a confirmation of COVID-19 infection.

PRINCIPLE OF THE TEST

The test consists of one test strip, which is integrated in a test cassette. This test strip consists of a highly specific neutralizing anti-SARS-CoV-2 antibody, coupled to colored particles (conjugate), and a membrane with one test line and one control line. The test line contains anti-SARS-CoV2 Spike glycoprotein (S1) monoclonal antibody, the control line consists of an antibody-binding protein. Test line and control line in the result window are not visible before applying a sample.

After the sample is pipetted into the sample well (S), the sample passes through the conjugate and the antigen in the samples bind to the conjugate. The antigen-conjugate complex migrates due to the capillary action to the site of the membrane where the monoclonal anti-SARS-CoV2 Spike glycoprotein (S1) antibody is immobilized (test line). If SARS-CoV-2 viral S1 antigens are present in the sample, they will bind to the test line and one colored line appears. The remaining complex migrates further across the membrane to the control zone ("C"). A colored line appears, indicating that the test was performed correctly.

SUPPLIED MATERIALS


Package sizes:

REF: COV19_AG_20 (20 Tests).

TEST COMPONENTS

- 20 Test cassettes, individually sealed in an aluminium bag with desiccant
- 20 empty specimen collection tubes
- 20 nozzle caps for the tubes
- 20 single use sterilized swabs for specimen collection, added according to Directive 93/42/EWG, manufacturer: Jiangsu Hanheng Medical Technology co., LTD., www.chinacytobrush.com Tel : +8613063969010; China, CE
- Extraction / dilution buffer in dropper vial, 14 mL
- 1 Tube stand (not shown)
- 1 Instructions for use

Note: Pictures may differ from the original.

MATERIALS NEEDED BUT NOT SUPPLIED

- Stopwatch

PREPARATION OF REAGENTS

All reagents are ready-to-use. No further preparation of reagents is necessary.

STABILITY AND STORAGE CONDITIONS

Store the test at 2 - 30°C. Unopened kit components (aluminium bags and buffer) are stable until the expiry date. The expiry date is printed on the labels of the aluminium bag, the buffer and the outer packaging. Do not use if the aluminum bag is damaged. **DO NOT FREEZE** or expose to temperatures above 30°C.

Aluminium pouch with test cassette:

Keep the test in unopened aluminium bag at 2 - 30°C!

Extraction / dilution buffer:

Store the buffer at 2 - 30°C. Unopened buffer is stable until the expiry date. After first opening the buffer is stable until the expiry date, if the bottle is tightly closed after every usage.

Single use sterilized swabs for specimen collection:

Store the swabs at 2 - 30°C. Do not use if the outer packaging is damaged!

WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use! Not for personal use!

Read the instructions carefully before performing the test.

- In accordance with Good Laboratory Practice (GLP), all laboratory devices employed should be regularly checked for the accuracy and precision.
- For professional in-vitro diagnostics only!
- Use all reagents within the expiry period (printed on the labels).
- Do not use reagents from different kit lots or batch codes and avoid mixing of reagents of different kit lots or batch codes.
- Avoid contamination of the reagents. Do not use the same container for several samples!
- Avoid repeated freezing and thawing of the samples because it could lead to denaturation of the antigens.
- Do not ingest or swallow! Do not eat, drink and smoke in the laboratory! Do not work without wearing protective clothing (gloves, safety glasses, safety mask and lab coat)! Avoid the contact of kit reagents with skin, eye or mucosa.
- All kit components should be considered as infectious agents. Decontaminate and dispose of residues of kit reagents and samples in accordance with local regulations, e.g. by autoclaving or using a disinfecting solution.
- Avoid touching of the membrane in the result window of the test device with your fingers (danger of contamination).
- Do not pipette samples and diluent directly onto the membrane in the result window of the test device.
- For single use only. The test is sensitive to moisture. Do not use if the outer packaging (aluminum bag) is damaged. After opening the aluminum bag, it must be used within 1 hour.
- All patients should be treated as potentially infectious. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens.
- Bring specimens to room temperature (preferably 15 - 30°C).
- Lipemic, hemolytic or bacterially contaminated samples must not be used.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.
- If infection with SARS-CoV-2 is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions and sent to state or local health departments for testing.
- Viral isolation in cell culture and initial characterization of viral agents recovered in cultures of SARS-CoV-2 specimens are NOT recommended, unless the works is done completely in a BSL3 laboratory using BSL3 practices and guidelines.

SAMPLE COLLECTION AND PREPARATION

Specimen Swab: (we recommend CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19²⁸).

Collecting from nasopharynx:

- Remove the swab from its packing.
- Insert flexible wire shaft swab through the nares parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient indicating contact with the nasopharynx.
- Gently, rub and roll the swab, rotate in a circular motion gently against the nasopharyngeal mucosa for 10 – 15 seconds (Rotating against the nasal wall).
- Leave the swab in place for several seconds to absorb secretions before removing. Withdraw the swab while making sure that the tip of the swab is wet (ensure the swab contains cells as well as mucus).

Collecting from oropharynx:

- Insert swab through the mouth into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, gums and remove the swab.
- Process the swab as soon as possible after collecting the specimen.

General Notes

- The test works best with fresh samples. Swabs specimens should be tested as soon as possible after collection. Use freshly collected specimens for best test performance.
- Avoid freezing and thawing specimens.
- If not tested immediately, swab specimens may be stored at 2-8°C for 24 hours after collection. Keep below -20°C for longer storage.
- Viability of some pathogens from specimens that were frozen and then thawed is greatly diminished and may result in false-negative test results.
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of sample with the interpretation of test results.


ATTENTION: Handle human nasopharynx and oropharynx as potentially infectious.

TEST PROCEDURE

Test cassette, buffer and patient's samples should be brought to room temperature (preferably 15 - 30°C) prior to testing. Do not open pouches until ready for testing.


Take the required number of test pouches sample collection tubes and swabs from the packaging kit. Open the aluminum pouch and place the cassette / s on a clean, non-absorbent flat surface. Label the test device with patient identification number. Place the sample collection tubes open (without nozzle cap) in the plastic rack.

After opening the aluminum bag, the test should be carried out within less than 1 hour as the test strip is sensitive to humidity.

1. Gently shake or vortex the extraction/dilution buffer bottle and drop 10 drops of the extraction/dilution buffer in the specimen collection tube.

2. Collect nasopharynx or oropharynx with the swab(s). **Attention! Use a new swab for each patient!** Follow the instructions in the section „Sample collection and preparation“ .

3. Insert the swab into the specimen collection tube. While squeezing the buffer tube, mix and stir the swab 10-15 times by compressing the walls of the tube against the swab. Keep the swap in the tube for 2 minutes.

4. Roll the swab head against the inner wall of the tube as you remove it. Try to release as much liquid as possible while squeezing the side of the tube and insert the nozzle cap into the sample specimen collection tube. Dispose of the used swab in compliance with your biohazard waste disposal procedure.

5. Carefully invert the extraction tube and hold the tube on the sample well "S". **Attention!** **Do not shake!**

6. Add 3 drops of extracted specimen into the sample well (S) by gently squeezing the tube. **Take care that the specimen is not dropped outside the sample well!**

7. Wait for 15 minutes and read the results directly from the test device (latest after 20 min). Do not read results after more than 20 minutes.

8. No test line indicates a negative result; two lines indicate a positive result.

INTERPRETATION OF RESULTS

- NEGATIVE:** Only one colored line appears in the control zone (control line "C", see Fig. 1A). In the test zone ("T") there should be no line visible.
- POSITIVE:** Two colored lines appear. One line should be visible in the control zone ("C") and one in the test zone ("T") (Fig. 1B and C). **The test lines "T" may be stronger or weaker than the control line "C".**
- DOUBTFUL:** Very weak shadow-like test line should be regarded as not clear. In this case it is recommended to take another sample from the same patient to measure it again using a new test.
- NOT VALID:** No control line visible and / or background color affects readability of test results.


Fig. 1: Schematic diagram of possible test results for **COVID-19 Ag Rapid Test**:

Negative result (A): only the control line appears; **Positive result (B) and (C):** two lines appear, test- and control line. **Invalid test (D):** only the test line appears.

QUALITY CONTROL

The **COVID-19 Ag Rapid Test** contains an internal control. A colored line in the control zone ("C") is considered as an internal procedural control. It confirms enough sample volume and correct test procedure. A clear background is an internal negative procedural control. If a background color appears in the result window and thereby the readability of the test results will be affected, the result may be invalid.

Insufficient sample volume or incorrect handling of the test are the most likely reasons for a lack of control line and / or a formation of background color which affects the readability of control / test lines. Check again the instructions of sample preparation and test procedure and repeat the test with a new test device. If the problem persists, contact the manufacturer or your local distributor.

PERFORMANCE CHARACTERISTICS

To estimate the reproducibility of the measurements, intra- and inter-assay variations as well as inter-operator variations and batch-to-batch variations were determined by measuring samples of different reactivity using **COVID-19 Ag Rapid Test**. No or insignificant intra- and inter-assay variations, inter-operator- and batch-to-batch variations were observed for **COVID-19 Ag Rapid Test**.

No **High dose hook effect** is observed for antigen concentrations up to 200000 ng/mL. Using samples with high SARS-CoV2 Spike S1 antigen concentration, it is observed that the control line becomes weaker. This was expected according to the test sandwich design, the control line was always visible on all tests. Additionally, results in all clinical positive samples also showed no High dose hook effect.

Diagnostic Sensitivity and Specificity:

For determination of clinical sensitivity / specificity, the results of **COVID-19 Ag Rapid Test** are compared by RT-PCR methods (both PCR comparators were performed according to the KIT manufacturer protocols. In the first PCR all reagents were from Thermo Fisher, the Platform is ABI Perkin Elmer and the detection carried out using real-time RT-PCR in accordance with the reference method published by Charité (Berlin)8. The second PCR comparator was according to the KIT manufacturer. The PCR Kit is the RealStar SARS-CoV-2 (altona diagnostics GmbH) and performed using the PCR cycler CFX 96 Touch (Bio-Rad Laboratories, Inc).

The results of clinical laboratory are taken as Gold standard (RT-PCR). All positive and negative samples were well defined and validated clinical samples provided by clinical labs. The results of PCR are taken as "true" (PCR positive or negative) for 242 samples (naso- and oropharyngeal swabs samples from symptomatic and asymptomatic patients). 126 samples (naso- and oropharyngeal swabs) with known of COVID-19 positivity and 116 negative samples (PCR negative swab samples, considered as healthy or other disease patients) were measured. This panel of samples with RT-PCR confirmed positive (Cycle threshold (Ct) range positive: 16.58 - 38.27) or negative patients was tested using **COVID-19 Ag Rapid Test**.

COVID-19 Ag Rapid Test showed **Diagnostic Sensitivity** of **94.44%** (CI = 88.89 to 97.74 %) at **Diagnostic specificity** of **98.28%** (CI = 93.28 to 99.79 %). The results of the comparative study are summarized in the following table (CI ="exact" Clopper-Pearson confidence intervals²⁰):

Group		Number of specimens	Clinical sensitivity in % (95% Confidence interval)	Clinical specificity in % (95% Confidence interval)
RT-PCR	positive (TP/FN)	126 (119/7)	94.44 (88.89 to 97.74)	98.28 (93.91 to 99.79)
	negative (TP/FN)	116 (114/2)		

The following table demonstrates the Correlation of the Diagnostic sensitivity results with different Ct values, which illustrate the clear dependency of the COVID-19 rapid test diagnostics sensitivity on the viral load and Ct values* of the used Gold standard/reference RT-PCR methods:

Ct Range	Clinical sensitivity in % (95% Confidence interval)
16.58 - 25	95.31 (86.91 to 99.02)
16.58 - 30	96.19 (90.53 to 98.95)
16.58 - 38.27	94.44 (88.89 to 97.74)
30 - 38.27	85.71 (63.66 to 96.95)

*Although the Ct values of PCR systems provide clues for the underlying virus concentration, please note that at the same virus concentration the Ct values may vary between different PCR methods. This depends on the sample volume used in the extraction, the proportion of the elution volume in the PCR approach and the extraction, elution and amplification efficiency differ between the different PCR methods.

Analytical sensitivity / Limit of detection (LoD):

The **Limit of detection (LoD)** was determined by COVID-19 Ag Rapid Test results out of the clinical investigations and with clinical samples panel with known Ct values. It was evaluated to be approximate at **Ct ≤ 30**.

COVID-19 Ag Test		
Samples expected results/ Ct values	% of positive results	% of positive results
Positive (Ct=27)	100.00	0.00
Positive (Ct=32)	66.66	33.33
Positive (Ct=23)	100.00	0.00
Positive (Ct=20)	100.00	0.00
Positive (Ct=31)	100.00	0.00
Positive (Ct=26)	100.00	0.00
Positive (Ct=23)	100.00	0.00
Positive (Ct=19)	100.00	0.00
Positive (Ct=19)	100.00	0.00
Positive (Ct=20)	100.00	0.00
Positive (Ct=25)	100.00	0.00
Positive (Ct=23)	100.00	0.00
Positive (Ct=26)	100.00	0.00
Positive (Ct=27)	66.66	33.33
Positive (Ct=22)	100.00	0.00
Positive (Ct=17)	100.00	0.00
Positive (Ct=30)	33.33	66.66
Positive (Ct=19)	100.00	0.00

Additionally, the Analytical sensitivity/Limit of detection (LoD) was determined by Rec. SARS-CoV2 Spike S1 spiked into a swab negative sample matrix with different concentration and has been found to be **≤ 0.1 ng/mL**.

Interfering Substances

Interference and cross-reactivity were assessed by respecting clinical data (from various similar test manuals and publications). To determine analytical specificity, samples were spiked by potential interfering substances.

Potentially interferents and cross-reacting substances/agents and are summarized below:

Neo-Synephrine (Phenylephrine)	10 %
Erythromycin	50 g/mL
Chloramphenicol	50 g/mL
Dopamine hydrochloride	10 g/mL
Human Albumin	110 mg/mL
Biotin	200 ng/mL
Caffeine	100 g/mL
Hemoglobin	1 mg/mL
Acetaminophen	50 g/mL
Acetylsalicylic acid	400 g/mL
Ibuprofen	400 g/mL
Zanamivir	5 mg/mL
Oseltamivirphosphat	10 mg/mL

No interference was observed for the substances tested.

Cross reactivity

Potentially cross-reacting organisms tested are summarized below (SARS-CoV-2 negative samples):

Virus/ Bakteriea/ Parasite (Strain)	Concentration
BCG ^{a)}	3×10^6 cells/mL
<i>Staphylococcus aureus</i> ^{a)}	3×10^6 cells/mL
<i>Mycobacterium smegmatis</i> mc ² 155 ^{a)}	3×10^6 cells/mL
<i>Mycobacterium avium</i> ^{a)}	3×10^6 cells/mL
Respiratory Viral Panel (Human Rhinovirus, Enterovirus und Adenovirus) ^{b)}	Positive swab sample
Influenza A ^{b)}	Positive swab sample
Influenza B ^{b)}	Positive swab sample
<i>Haemophilus influenzae</i> ^{b)}	Positive Sputum sample
<i>Streptococcus pneumoniae</i> ^{b)}	Positive Sputum sample
Yeast ^{b)}	Positive Sputum sample
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) ^{b)}	Positive Sputum sample
<i>Pseudomonas aeruginosa</i> PAO1 ^{a)}	3×10^6 cells/mL
<i>E. coli</i> ^{a)}	3×10^6 cells/mL
<i>Mycobacterium paratuberculosis</i> 032045 ^{a)}	3×10^6 cells/mL
<i>Mycobacterium marinum</i> 3-1521/68 ^{a)}	3×10^6 cells/mL
<i>Mycobacterium kansassii</i> ^{a)}	3×10^6 cells/mL
<i>Mycobacterium gordoneae</i> ^{a)}	3×10^6 cells/mL
<i>Mycobacterium vaccae</i> NC 10916 ^{a)}	3×10^6 cells/mL
Human Coronavirus 229E (HCoV-229E) ^{c)}	$3,30 \times 10^6$ TCID ₅₀ /mL
Human Coronavirus NL63 (HCoV-NL63) ^{c)}	$1,00 \times 10^4$ TCID ₅₀ /mL
Human Coronavirus OC43 (HCoV-OC43) ^{c)}	$3,70 \times 10^2$ TCID ₅₀ /mL
MERS-CoV Spike S1 Protein in Negativ-Probenmatrix ^{d)}	1 g/mL
	500 ng/mL
	250 ng/mL

Description (SARS-CoV-2 negative samples):

- a) Cell culture spiked in swab Negative sample matrix.
- b) Sample matrix from Discovery Life Sciences, AL, 35806 USA.
- c) Inactivated virus culture from National Infection Service, Public Health England (PHE).
- d) Recombinant protein from Sino Biological Europe GmbH

There was no cross-reaction with potential cross-reactive organisms/substances tested as shown above.

In addition, the monoclonal detection antibody used has shown no cross-reactivity in ELISA with SARS-CoV Spike S1-mFc Protein, SARS-CoV Spike RBD-His Protein, MERS-CoV Spike S1 Protein, HCoV-HKU1 (isolate N1) Spike S1 Protein, HCoV-HKU1 (isolate N5) Spike S1 Protein, HCoV-NL63 Spike S1 Protein, HCoV-229E Spike S1 Protein, HCoV-OC43 Spike S1+S2 ECD Protein (according to the manufacturers specifications).

Inclusivity (analytical sensitivity):

The target of **COVID-19 Ag Rapid Test** is SARS-CoV-2 (2019-nCoV) surface glycoprotein (Spike protein-S1). In silico analysis of published different SARS-CoV-2 strain sequences showed that 100% of analyzed published SARS-CoV-2 strains will be detectable using this assay target protein and our novel neutralizing antibody.

As shown in table below the analysis showed very high homology of more than 99.5%.

Strain	NCBI GenBank ACCESSION Nr.	% homology
Wuhan-Hu-1	MN908947.3	100%
SARS-CoV-2/Hu/Kng/19-437	LC534419.1	100%
SARS-CoV-2/Hu/DP/Kng/19-031	LC534418.1	100%
SARS-CoV-2/Hu/DP/Kng/19-020	LC528232.1	100%
FDAARGOS_983(isolate="SARS-CoV-2/human/USA/USA-WA1/2020)	MT246667.1	100%
SARS-CoV-2/human/PAK/KP-RMI-01/2020	MW242667.1	100%
SARS-CoV-2/human/POL/PL_MCB_10/2020	MW273792.1	100%
SARS-CoV-2/human/JPN/UT-NCGM02/2020	MW219695.1	100%
SARS-CoV-2/human/KHM/Kunming_kms-2/2020	MW341443.1	100%
SARS-CoV-2/human/RUS/1150/2020	MW332225.1	99.85%
SARS-CoV-2/human/BHR/341036861/2020	MW345922.1	99.85%
SARS-CoV-2/human/IND/GBRC417b/2020	MW242689.1	99.85%
SARS-CoV-2/human/AUS/VIC17053/2020	MW320822.1	99.70%
SARS-CoV-2/human/SRB/KV0052-12-05/202	MW266938.1	99.70%
SARS-CoV-2/human/USA/GA-CDC-8142/2020	MW343786.1	99.70%
SARS-CoV-2/human/USA/GA-CDC-7877/2020	MW343784.1	99.70%

LIMITATIONS

Follow the instructions of the test procedure and interpretation of results carefully!

- Follow the instructions of the test procedure and interpretation of results carefully! Insufficient sample volume or incorrect handling of the test procedure are the most likely reasons for not reaching the required QC criteria of test performance (see section “Quality control of test”).
- The **COVID-19 Ag Rapid Test** intended for testing a swab directly without elution in viral transport media as dilution will result in decreased detection of low positive samples that are near the limit of detection of the test. **Therefore, Swab samples eluted in VTM are not appropriate for use in this test.**
- A NEGATIVE result does not preclude SARS-CoV-2 infection and should be confirmed via molecular assay. Note that questionable results require further confirmation. If the result is not clear, another sample should be taken from the same patient and checked again.
- **COVID-19 Ag Rapid Test** is for professional in vitro diagnostic use and should only be used for the qualitative detection of SARS-CoV-2 antigen. The intensity of color in a positive band should not be evaluated as “quantitative or semi-quantitative”.
- **Both viable and nonviable SARS-CoV-2 viruses are detectable with the COVID-19 Ag Rapid Test.** As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- The performance of the **COVID-19 Ag Rapid Test** was evaluated using the procedures provided in this product insert only.
- Modifications to these procedures may alter the performance of the test.
- Negative results should be treated as presumptive and tested with an alternative authorized molecular assay, if necessary, for clinical management, including infection control.
- False negative results may occur if a specimen is improperly collected, transported, or handled. Negative results should be considered in the context of a patient’s recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

-
- As with any molecular test, mutations within the antigen target regions of the **COVID-19 Ag Rapid Test** could affect the binding resulting in failure to detect the presence of the virus.
 - The **COVID-19 Ag Rapid Test** does not replace the direct detection by PCR. It is important to note that a positive **COVID-19 Ag Rapid Test** result against SARS-CoV-2 indicates that an infection has taken place.
 - It is recommended to consider the results of the test in combination with the clinical status of each patient, the results of other diagnostic tests and the epidemiological background information. If a patient sample has tested positive, further confirmatory tests should be performed (e.g. PCR, clinical symptoms). For a final diagnosis, include all available information on a given patient.
 - The **COVID-19 Ag Rapid Test** does not assess the immune response and for this you need other testing methods **available separately from LIONEX** like **COVID-19 IgG/IgM Rapid Test** or **COVID-19 ELISA (Human IgG) / COVELISA® COVID-19 IgM ELISA**.
 - As with all diagnostic tests, a definitive clinical diagnosis should not be based solely on the results of a test but should only be made by the physician based on an evaluation of clinical and laboratory findings.
 - **Although the Ct values of PCR systems provide clues for the underlying virus concentration, they can only be compared with restrictions between different PCR methods, because the sample volume used in the extraction, the proportion of the elution volume in the PCR approach and the extraction, elution and amplification efficiency differ between the different PCR methods.**

INSTRUCTIONS FOR USE

Disposable sampling swab


【Product Name】 Disposable sampling swab

【Product Model】 Type A(A-01, A-02, A-03, A-04, A-05), Type B, Type C, Type D, Type E, Type F

【Product structure】

1. The Type A sampling swab consists of a sampling head and a sampling handle. The sampling head is a fluffy head made of medical grade PA66 (polyamide 66) material. The sample handle material is made of ABS (Acrylonitrile Butadiene Styrene) or PP (Polypropylene) material. Type B and Type C consist of a sampling swab, sheath, and protective cap; the sheath and cap are made of PP (Polypropylene) material. The type D, E, F sampling swab consists of a sampling head and a sampling handle. The sampling head of type D is made of Polyester, The sampling head of type E is made of Rayon, The sampling head of type F is made of Absorbent cotton. The sampling handle of type D, E, F is made of PS(Polystyrene) or PP or wood.

2. The structural schematic of the sampling swab are shown in following Figures.

【Intended Use】

It is mainly used for clinical cell sampling and smear.

【Contraindication】

No

【Steps to use】

1. Open the package, gently insert the sample swab into the sampling site while holding the handle.
2. Gently rotate the sampling swab 3~5 turns and slowly remove it.
3. The extracted sample is placed in the sample collection tube, and then sealed to complete the sampling. The B, C-type can be directly covered with a sheath after sampling.

【Precaution】

1. Please check the packaging and the product carefully before use. If either of them is damaged, it is strictly forbidden to use.
2. This product is for one-time use only. It should be strictly in accordance with the aseptic operation specifications. Please destroy it immediately after use. It is forbidden to reuse. The disposal of wastes shall be carried out in accordance with national environmental protection laws and regulations.
3. Please read the instruction manual carefully, pay attention to the expiry date of the product before use. Do not use it after the deadline.

【Definitions of Signs】

	Sterilized using ethylene oxide		Do not re-use
	Do not use if package is damaged		Do not re-sterilize
	Caution		Fragile, handle with care
	Upward		Keep dry
	Batch code		Use-by date

- 3 -

	Manufacturer		Authorized representative in the European Community
	Date of manufacture		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Stacking limit by number		"CE" Marking

【Storage and transportation】

1. The product should be stored in a room with a temperature of -10°C~40°C, relative humidity of not more than 80%, good ventilation, no corrosive gas.
2. The product should be protected from heavy pressure, strong vibration, direct sunlight, high temperature, moisture, rain and corrosive gases or liquids during transportation.

【Sterilization】 EO sterilization

【Shelf Life】 5 years after sterilization

Jiangsu HanHeng Medical Technology Co., Ltd.
16-B4, #1 North Qingyang Road, Tianning District,
213017 Changzhou, Jiangsu, China

Name: Luxus Lebenswelt GmbH
Add: Kochstr. 1, 47877, Willich, Germany
Contact Person: Liu Sun
Tel/Fax: 0049-1715605732
E-mail: Info.m@luxuslw.de

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【Sterilisation】

【Lagerfähigkeit】

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Jiangsu HanHeng Medical Technology Co., Ltd.
16-B4, #1 North Qingyang Road, Tianning District,
213017 Changzhou, Jiangsu, China
Name: Luxus Lebenswelt GmbH

Adresse: Kochstr. 1, 47877, Willich, Germany
Kontaktperson: Liu Sun
Tel/Fax: 0049-1715605732
E-mail: Info.m@luxuslw.de

【Version】 2020.09.10

A/2

LITERATURE / LITERATUR

- [1] Weiss SR, Leibowitz JL. Novel coronavirus (2019-nCoV), World health Organisation (WHO), 2020. Coronavirus pathogenesis. *Adv Virus Res* 2011;81:85-164. PMID:22094080 DOI:10.1016/B978-0-12-385885-6.00009-2
- [2] World Health Organization (WHO). WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. [Accessed 26 Jan 2020]. <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>
- [3] World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>
- [4] WHO (2020). Novel Coronavirus (2019-nCoV) Situation Report 23. https://www.who.int/docs/default-source/coronaviruse/situation-reports/20200212-sitrep-23-ncov.pdf?sfvrsn=41e9fb78_4.
- [5] Zhou, P., Yang, X.L., Wang, X.G., Hu, B., Zhang, L., Zhang, W., Si, H.R., Zhu, Y., Li, B., Huang, C.L., et al. (2020). A pneumonia outbreak associated with anew coronavirus of probable bat origin. *Nature*. <https://doi.org/10.1038/s41586-020-2012-7>.
- [6] Menachery, V.D., Dinnon, K.H., III, Yount, B.L., Jr., McAnarney, E.T., Gralinski, L.E., Hale, A., Graham, R.L., Scobey, T., Anthony, S.J., Wang, L., et al. (2020). Trypsin treatment unlocks barrier for zoonotic bat coronaviruses infection. *J. Virol.* 94 <https://doi.org/10.1128/JVI.01774-19>.
- [7] Iwata-Yoshikawa, N., Okamura, T., Shimizu, Y., Hasegawa, H., Takeda, M., and Nagata, N. (2019). TMPRSS2 Contributes to Virus Spread and Immunopathology in the Airways of Murine Models after Coronavirus Infection. *J. Virol.* 93 <https://doi.org/10.1128/JVI.01815-18>.
- [8] Corman, V.M., Landt, O., Kaiser, M., Molenkamp, R., Meijer, A., Chu, D.K., Bleicker, T., Bru" nink, S., Schneider, J., Schmidt, M.L., et al. (2020). Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *Euro Surveill.* 25 <https://doi.org/10.2807/1560-7917.ES.2020.25.3.2000045>.
- [9] Hoffmann et al., SARS-CoV-2 Cell Entry Depends on ACE2 and TMPRSS2 and Is Blocked by a ClinicallyProven Protease Inhibitor, *Cell* (2020), <https://doi.org/10.1016/j.cell.2020.02.052>.
- [10] Daniel Wrapp, Nianshuang Wang, Kizzmekia S. Corbett, Jory A. Goldsmith, Ching-Lin Hsieh, Olubukola Abiona, Barney S. Graham and Jason S. McLellan. Cryo-EM structure of the 2019-nCoV spike in the prefusion conformation. originally published online February 19, 2020 DOI: 10.1126/science.abb2507 (6483), 1260-1263. 367 *Science*
- [11] Corman, V.M., Landt, O., Kaiser, M., Molenkamp, R., Meijer, A., Chu, D.K., Bleicker, T., Bru" nink, S., Schneider, J., Schmidt, M.L., et al. (2020). Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *Euro Surveill.* 25 <https://doi.org/10.2807/1560-7917.ES.2020.25.3.2000045>.
- [12] Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019;17:181-192. PMID:30531947 DOI:10.1038/s41579-018-0118-9
- [13] Gallagher and Buchmeier (2001). Coronavirus spike proteins in viral entry and pathogenesis. *Virology*. 279(2):371-4.
- [14] Ji et al. (2020). Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human. *J Med Virol.* 2020;10:1002/jmv.25682. doi:10.1002/jmv.25682.
- [15] Li F. (2016). Structure, Function, and Evolution of Coronavirus Spike Proteins. *Annu Rev Virol.* 3(1):237-261.
- [16] Lu et al. (2015). Bat-to-human: spike features determining 'host jump' of coronaviruses SARS-CoV, MERS-CoV, and beyond. *Trends Microbiol.* 23(8):468-78.
- [17] Lu R, Zhao X, Li J, et al. (2020). Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *Lancet.* S0140-6736(20)30251-8. doi:10.1016/S0140-6736(20)30251-8.
- [18] Song et al. (2018). Cryo-EM structure of the SARS coronavirus spike glycoprotein in complex with its host cell receptor ACE2. *PLoS Pathog.* 2018 Aug; 14(8).
- [19] Su et al. (2016). Epidemiology, Genetic Recombination, and Pathogenesis of Coronaviruses. *Trends Microbiol.* 2016 Jun; 24(6):490-502.
- [20] Clopper, C.; Pearson, E. S. (1934). "The use of confidence or fiducial limits illustrated in the case of the binomial". *Biometrika*. 26: 404–413. doi:10.1093/biomet/26.4.404.
- [21] Lee CY-P, Lin RTP, Renia L and Ng LFP (2020) Serological Approaches for COVID-19: Epidemiologic Perspective on Surveillance and Control. *Front. Immunol.* 11:879. doi: 10.3389/fimmu.2020.00879.
- [22] Long QX, Deng HJ, Chen J, Hu J, Liu BZ, Liao P, et al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. *medRxiv [Preprint]*. (2020). doi: 10.1101/2020.03.18.20038018
- [23] Lee N, Chan PK, Ip M, Wong E, Ho J, Ho C, et al. Anti-SARS-CoV IgG response in relation to disease severity of severe acute respiratory syndrome. *J Clin Virol.* (2006) 35:179–84. doi: 10.1016/j.jcv.2005.07.005.
- [24] DC Emergency Operations Center (EOC), National Center for Immunization and Respiratory Diseases (NCIRD), Division of Viral Diseases. Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19. Last Updated July 8, 2020. <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>.
- [25] Wyllie AL, Fournier J, Casanovas-Massana A, Campbell M, Tokuyama M, Vijayakumar P, et al. Saliva or Nasopharyngeal Swab Specimens for Detection of SARS-CoV-2. *New England Journal of Medicine* (2020) 383(13):1283-6. doi: 10.1056/NEJMc2016359.
- [26] Tang Y-W, Schmitz JE, Persing DH, Stratton CW. 2020. Laboratory diagnosis of COVID-19: current issues and challenges. *J Clin Microbiol* 58:e00512-20. <https://doi.org/10.1128/JCM.00512-20>.
- [27] Mak GC, Cheng PK, Lau SS, Wong KK, Lau CS, Lam ET, et al. Evaluation of rapid antigen test for detection of SARS-CoV-2 virus. *J Clin Virol.* 2020;129:104500.
- [28] Bruning AHL, Leeflang MMG, Vos J, Spijker R, de Jong MD, Wolthers KC, et al. Rapid Tests for Influenza, Respiratory Syncytial Virus, and Other Respiratory Viruses: A Systematic Review and Meta-analysis. *Clin Infect Dis.* 2017;65(6):1026-32.
- [29] Mertens P, De Vos N, Martiny D, Jassoy C, Mirazimi A, Cuypers L, et al. Development and Potential Usefulness of the COVID-19 Ag Respiratory Strip Diagnostic Assay in a Pandemic Context. *Front Med (Lausanne)*. 2020;7:225.
- [30] Dinnes J, Deeks JJ, Adriano A, Berhane S, Davenport C, Dittrich S, et al. Rapid, point-of-care antigen and molecular-based tests for diagnosis of SARS-CoV-2 infection. *Cochrane Database Syst Rev.* 2020;8:CD013705.
- [31] Krüger LJ, Gaeddert M, Köppel L, Brümmer LE, Gottschalk C, Miranda IB, et al. Evaluation of the accuracy, ease of use and limit of detection of novel, rapid, antigen-detecting point-of-care diagnostics for SARS-CoV-2. *medRxiv*. 2020:2020.10.01.20203836.

- [32] Ghebremedhin B, Engelmann I, Konig W, Konig B. Comparison of the performance of the rapid antigen detection actim Influenza A&B test and RT-PCR in different respiratory specimens. *J Med Microbiol.* 2009;58(Pt 3):365-70.
- [33] Lindner AK, Nikolai O, Kausch F, Wintel M, Hommes F, Gertler M, et al. Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab. *medRxiv* (2020).
- [34] Corman VM, Haage VC, Bleicker T, Schmidt ML, Mühlmann B, Zuchowski M, et al. Comparison of seven commercial SARS-CoV-2 rapid Point-of-Care Antigen tests. *medRxiv* (2020).

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CONTACT INFORMATION / KONTAKTINFORMATIONEN

For more information and technical assistance, please contact /

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